

# Datasheet ARKY indoor cabinet



ARKY cabinets are designed and manufactured by Defibrion BV in Groningen, the Netherlands. The ARKY indoor cabinet (article 60.111, 60.112, 60.121, 60.122 and 60.123) is suitable for the indoor storage of an Automated External Defibrillator (AED).

## SPECIFICATIONS

### Size and weight

|            |                   |
|------------|-------------------|
| Dimensions | 380 x 380 x 200mm |
| Weight     | 5kg               |

### Electronic features

|              |                               |
|--------------|-------------------------------|
| Alarm        | Acoustic                      |
| Power supply | 3x Button Cell Batteries LR44 |
| Output       | 105 dB                        |

### Materials used

|               |                |
|---------------|----------------|
| Cabinet shell | Metal          |
| Paint         | Powder coating |
| Hinge         | Aluminium      |
| Window        | Plastic        |

### Operating environment

|                                 |                  |
|---------------------------------|------------------|
| Operating / storage temperature | From 5°C to 30°C |
|---------------------------------|------------------|

### Options

|        |                |
|--------|----------------|
| Colors | Green or white |
|--------|----------------|