

# Martek CardioSAFE™

## Your Livesaving Coach

### Always Ready

Without any need for electricity or batteries, the Martek CardioSAFE™ is entirely mechanical – meaning it's always ready to save a life.

### Simple Instructions

Easy-to-read instructions are printed on the Martek CardioSAFE™, so any bystander can use it with confidence in the case of a cardiac emergency.

### Perfect Positioning

The Martek CardioSAFE™ is designed to sit perfectly in the correct positioning on a patient's chest, so cardiac massage is always delivered exactly where it's needed – and at the right depth for true effectiveness.

### Audio Feedback

A special spring creates an audible feedback signal when the perfect amount of pressure is applied to the chest. It will also trigger when the spring is released, enabling the rescuer to reach the required frequency of approximately 120 compressions per minute for successful CPR.



Data		Body Materials	
Length	16.6 cm	LEXAN Polycarbonate [CASRN 25971-63-5]	
Width	11.4 cm	Polyurethan flex foam KO T9020H	
Height	5.8 cm		
Weight	125 g		
Pressure		Certificates	
High end stainless steel springs		CE certified	
410 N/m		FDA 510 (k) pending	

**martek**™  
LIFECARE  
Making Life Better