









Tough Polyethylene body

- Corrosion resistant
- > Impact resistant
- Good thermal insulation
- Lightweight

Polycarbonate Door

- > Same material used to make riot shields
- > Tough & lightweight

Metal hinges with security fixings

> Minimise the chances of damage or vandalism

Double Skin Construction

> Tough and thermally efficient

Low Energy Consumption

Inexpensive to run

3 Bolt fixing to wall

> Ensures a stable fixing even on uneven surfaces

Protection against the elements

- > IP66 rated against dust and water ingress
- Thermostatically controlled heater prevents freezing
- > Air circulation fan prevents condensation
- > Safety cut out

LED Lighting

- Efficient internal illumination
- > Highly visible at night time

Unique Location code

Allows Emergency Services to pin-point your location



Stainless Steel Mechanical Lock

- > High quality, corrosion resistant
- In-built clutch to prevent forced entry damage to lock
- > Two year manufacturer's warranty
- Single handed entry
- > Large handle easier to open

Option for bespoke company logo on front of cabinet at an additional charge

Lockable and non-lockable versions available

For more information visit: defibsafe.co.uk Email: info@defibsafe.co.uk Call: 01252 344007













DefibSafe2 Cabinet: Lockable and Unlocked

Product Codes

Product Code: 91001 (Lockable) **Product Code:** 91006 (Unlocked)

Dimensions

APPLIES TO BOTH CABINETS

External Dimensions

Height: 625mm

Width: 440mm

Depth: 245mm

Depth including Lock Handle: 285mm

Internal Dimensions

Floor to shelf – 330mm

Width - 285mm

Depth with door closed – 160mm

Product Package

Cabinet; Installation Manual; Installation Template; 3 x washers; 1 x PVC Conduit Pipe; 3 x Expansion Bolts; 1 x UK Plug



Lockable (with code pad)



Unlocked



Interior

For more information visit: defibsafe.co.uk Email: info@defibsafe.co.uk Call: 01252 344007